



IMSF Workgroup Seminar Series - Instrumentation

Harnessing the power of the unreduced data in FTMS or how to squeeze out maximum information from ion signals

وداعدته بالقال إدبريا أعالي وروالا أنوروس أيروا فالبامير راونان معاديا روافيان وراك المراز والكريمونان

16 March 2021 @ 17:00 CET

Spectroswiss Lausanne, Switzerland

Yury Tsybin tsybin@spectroswiss.ch

FTMS: Fourier Transform Mass Spectrometry

- Ion identity (m/z) is encoded as a frequency of ion oscillations in an ion trap ullet
- Frequencies of ion oscillations are measured as time-domain signals (transients) ullet
- Fourier transform (FT) decodes transients to reveal frequencies (m/z) values ullet



Orbitrap[™] families: LTQ Orbitrap; Exactive; Q Exactive; Exploris; Fusion; Lumos; Eclipse

FTMS is the highest resolution & mass accuracy MS technique



Magnetic field-based Ion Cyclotron Resonance (ICR) Magnetic Resonance MS (MRMS)





The use of the unreduced data aims to provide performance improvement



metabolites identified & quantified

proteins identified & quantified

% sequence coverage in protein structural analysis

confidence for precursor and product ion detection

The Unreduced Data in FTMS

- I. The Unreduced Data in FTMS
- II. Data Reduction Approaches
- III. How to Get and Process the Unreduced Data
- IV. Examples of Applications



The Unreduced Data in FTMS

I. The Unreduced Data in FTMS

- II. Data Reduction Approaches
- III. How to Get and Process the Unreduced Data
- IV. Examples of Applications





- $R \sim T$ Resolution increases linearly with transient length (detection period) $S/N \sim \sqrt{T}$
 - Sensitivity increases as a sqrt of transient length (detection period)

Transients & aFT mass spectra provide equal information output



Time-Domain Ion Signal: Transient The Unreduced Data



example: two waves of the same frequency and amplitude, but different phase





Absorption Mode FT (aFT) Mass Spectrum The Unreduced Data



Simulations: IMSF, [M+H]⁺, 497.2 *m/z*, Orbitrap Q Exactive, 15k @ 200 *m/z*



aFT preserves all information:

- Frequency
- Amplitude
- Phase

noise (+ and -)



Simulation of FTMS data: DOI: 10.1021/jasms.0c00190

Why Reduce FTMS Data?

- The unreduced data file size:
 - a single transient: 1-20 MB;
 - a single unreduced mass spectrum: 2-40 MB
 - a single imaging/LC-MS run: **1-50 GB**
- Technological challenges/cost
 - High-performance electronics for transient acquisition
 - Real-time digital signal processing algorithms (on a chip)
 - Efficient software tools for post-processing of "big data"
- Ease of use
 - Extended data processing time for the unreduced data processing
 - A larger storage space is needed for the unreduced data





hip) ta"

ta processing data

9

Cost of Reducing FTMS Data?

- It is a one way street lacksquare
 - Once data is reduced, the full data cannot be recovered
 - (Unless the original data is also saved defeating the purpose)
- Data reduction always loses information/capability
 - Spectral quality •
 - Unresolved components
 - Noise distribution
 - Unexpected analytes
 - Sub-threshold peaks
 - Etc.

Note: the reduced data use is sufficient for many FTMS applications









The Unreduced Data in FTMS

I. The Unreduced Data in FTMS

II. Data Reduction Approaches

- III. How to Get and Process the Unreduced Data
- IV. Examples of Applications



Data Reduction Approaches in FTMS

- Transient length: ion detection during a reduced period of time \bullet (*not* all time ions induce signals)
- Calibration of initial phases: post-processing instead of in real time
- Performing other modes of FT: magnitude (mFT) or enhanced (eFT)
 - Noise reduction/thresholding: full vs. reduced profile
 - Visible mass (m/z) range: low and high mass cut-off
 - Centroiding mass spectra
 - Deisotoping, deconvolution
 - Transients are not stored, typically







Tsybin et al., DOI: 10.1016/B978-0-12-814013-0.00005-3



The Unreduced Data in FTMS

- I. The Unreduced Data in FTMS
- II. Data Reduction Approaches

III. How to Get and Process the Unreduced Data

IV. Examples of Applications



How to Get & Process the Unreduced FTMS Data







ExplorisTM schematics: https://planetorbitrap.com/

Overcoming the technological limitations in acquiring phased & full transients



JASMS 2020, 31, 257-266

DOI: 10.1021/jasms.9b000

FTMS Booster: a high-performance DAQ system

- Analogue transient capture after the original built-in pre-amplifier
- Phase artefact-free digitization (sampling) of transients
- Ability to maximize the duty cycle of ion detection event (throughput)
- Ability to record (much) longer transients (due to enhanced flexibility)
- Transient storage (.H5 files) for post-processing, reduced file size



FPGA: field programmable gate array (chip)



e to enhanced flexibility) ng, reduced file size

JASMS 2020, 31, 257-266



m/z



Spectroswiss Software Requirements for Processing of the «Big Data»

| Netu Navigate Help Import Data | II: TTR: SRUUK JIII/ARV IZ. GRUpHe new Re(DK) | | | naskavnask | |
|-----------------------------------|--|--|--|--|--|
| Data Overview Data Processing | FTMS Isotopic Simulator | | | ETMS DATA ANALYSIS | |
| Tools | FTMI isotopi smuldior: 🖌 Epuni kovi – Hole polic 🖌 Hole polic 🖌 Hole balles | | | | |
| | Shter de | | | | |
| | Bobot retings Expans Ners, Rangate Hor | nies suz | | | |
| | V Ferrivis (sectral): cz Import Data Data Overview | 19. Generate Mi Intege | | peak <mark>by</mark> peak | |
| | Data Processing Machiner(i) Post-processing | Generate MS image | | FTMS DATA ANALYSIS | |
| | e in a surges in Visualize Results Tacis Innection make: ja Denry date: nu Denry | Plot an image at specified mase-to-charge ratio(s) for an imput data file from | DESI/MALDI MSI experiment. | | |
| | | | | | |
| | | Constant Terrate Contractor | | | |
| | | Tel mass of a wedge ensember 8 Bit image of m Right according to | | | |
| | | 1K: 104 41 106 47 106 43 106 44 106 45 106 45 106 44 10 | 41 100-45 100-471 100-471 11 14 15 16 | 17 11 14 20 sil mages | |
| | | Libite Denity | | | |
| | | | | | |
| | | nghe 10.5410000 mg/s 200.0000 mg/s 20.0000 mg/s 20.00000 mg/s 20.0000 mg/s 2 | | | |
| | forepound litt. colort we from (| A WELSE | | | |
| | Lonous Relaites Fore | CONTRACTOR OF A CONTRACTOR | ANTI STREAME AND | and the second sec | |
| | 1.4 C +681 G254 M | average and a second | NUMBER STREET | - LANDOURS | |
| | La Constantia | And A Constanting BASSARA | | | |
| | 1.1 C 04081 (255 M | | | | |
| | 1.4 C (HIBL CZ3) H | a start range and a second | | Juran Shifteraye | |
| | 1.4 \$ -9665 (259-94 1.6 \$ 0608 (254-94 | H111 122 40050 | mate: 702.9633(78H1) | The State St | |
| | 2.5 C +HONO C1914 | | | | |
| | 773 \$ AUD45 C7414 | | | | |
| | | AND RES STATEMENT AND RESSON | | 12: E4 66/1077 | |
| | | | | | |
| | | | | | |
| | | | Image generation: 300.000% (elapsed time: 0:00:22 | 0 | |
| | | | | | |
| Si | mulated | l transient | FT mas | s spectrum | |
| | | | | | |
| | 1- | | | \bigwedge | |
| | | | | | |
| | Accessed to constitute of | There | ET | | |
| | | | FI | | |
| | | | | | |
| | | T | | | |
| 100 | | | | | |
| - | | | | | |
| THE | | = | | | |
| : | | | | ana Trana Lana Tana Pana . | |
| | | | | | |
| 1 | ··· I | | | | |
| | | | | | |
| | MHz | Time | | m/z | |

- Processing of any size and number ${\bullet}$ transients: user-defined FT parameters
- Mass spectra generation in absorption mode ulletFT (aFT): vs. the eFT & mFT (e.g., RAW)
- Processing of any size & number of mass spectra (*e.g.*, RAW in reduced & full profile)
- Processing of mass spectra (*e.g.*, RAW) from \bullet different experiments (*e.g.*, tech. replicates)
- Transient-mediated accurate simulation FTMS data (*DOI: 10.1021/jasms.0c00190*)
- Data analysis of mass spectra (*e.g.*, RAW) lacksquareusing the accurately simulated FTMS data

Peak-by-Peak software package (Spectroswiss) fulfills these requirements

Of

of

Another Example of the «Big Data» FTMS Software



- Suite of MS data (post) processing modules •
- Uses transient as starting point unreduced data! ullet
- Can be combined into many different workflows
 - AutoPhaser aFT mode processing
 - AutoPiquer peak picking
 - AutoSeequer top-down assignment
 - **Discharger** deconvolution & decomposition
 - AutoLogis DOM and petroleomics:



Streamlined recalibration



aFT image reprocessing



TOF import

aFT fine structure

Deconvolution Decomposition

aFT peak detection

https://www.kilgourlab.com/





David Kilgour

10.1016/j.marchem.2021.103955





The Unreduced Data in FTMS

- I. The Unreduced Data in FTMS
- II. Data Reduction Approaches
- III. How to Get and Process the Unreduced Data

IV. Examples of Applications



Examples of Applications

Small molecules (metabolites/peptides)

Data averaging: sensitivity & quantitation accuracy

Full transients: ultra-high resolution, throughput

Transient post-processing: super-resolution algorithms

single ion counting





Large molecules (proteins)

Data averaging: top-down & middle-down sequencing

Transient post-processing: enhanced spectral quality

Transient post-processing:

Examples of Applications

Small molecules (metabolites/peptides)

Data averaging: sensitivity & quantitation accuracy

Full transients: ultra-high resolution, throughput

Transient post-processing: super-resolution algorithms

single ion counting





Large molecules (proteins)

Data averaging: top-down & middle-down sequencing

Transient post-processing: enhanced spectral quality

Transient post-processing:



Srzentic et al., Anal. Chem. (2018) 12527

Enhanced Protein Sequencing with Top/Middle-Down 2D Data Averaging: Use of Multiple LC-MS/MS Technical Replicates



Srzentic et al., Anal. Chem. (2018) 12527; Fornelli et al., J. Proteomics (2017) 67



24

Examples of Applications

Small molecules (metabolites/peptides)

Data averaging: sensitivity & quantitation accuracy

Full transients: ultra-high resolution, throughput

Transient post-processing: super-resolution algorithms

single ion counting





Large molecules (proteins)

Data averaging: top-down & middle-down sequencing

Transient post-processing: enhanced spectral quality

Transient post-processing:

Data Post-Processing for Improved Protein Analysis Simulations, Orbitrap D20, 5 kV Ovalbumin (~43 kDa), protonated, z=38 $60k@200 m/z, T_{acq} = 128 ms$ ¹³C₂₆¹⁵N₂¹⁸O₂³⁴S₂ ¹³C_o¹⁵N Ovalbumin, z=38, including noise (gaussian) а b Amplitude, a.u. (x1e+12) T_1 T_2 T_4 T_3 13C815N2 Intensity, a.u. (x1e+09) 10 20-1 1.1 li itt т L 1128.306 1128.311 1128.996 1128.997 0-10 4+17 A+31 -10 0.02 0.04 0.1 0.12 0.06 0.08 n 1128.8 1129 1128.4 1128.6 Time, s m/z T_3 300 12 С d - mFT, full window Signal-to-noise ratio (SNR) Intensity, a.u. (x1e+12) = 8 ms 8 $T_1 = 8 \text{ ms}$ 250 aFT, half window 8 - aFT, full window $T_2 = 16 \text{ ms}$ $T_2 = 16 \text{ ms}$ 6 200 $T_3 = 32 \text{ ms}$ $T_3 = 32 \text{ ms}$ 6 4 150 $T_4 = 64 \text{ ms}$ $T_4 = 64 \text{ ms}$ 2 100 0 50 -2 0 0.02 0.06 0.08 0.1 0.04 0 1132 1133 1133 1132 m/z m/zDetection period, s

Experiment: Fornelli et al., J. Proteome Res. 2017, 16, 609–618





Examples of Applications

PULL PEAKS

NPAR7

Small molecules (metabolites/peptides)

Data averaging: sensitivity & quantitation accuracy

Full transients: ultra-high resolution, throughput

Transient post-processing: super-resolution algorithms

Large molecules (proteins)

single ion counting



- Data averaging: top-down & middle-down sequencing
- Transient post-processing: enhanced spectral quality
- Transient post-processing:

Single Ion Counting: Resolving Protein Interference



STORI: Selective Temporal Overview of Resonant Ions, DOI: 10.1021/jasms.8b06253



Detection of single ions reduces space artefacts charge (peak interference), hides but peak's charge state

Charge state of a single ion peak can be deduced from its transient, *e.g.*, via STORI (I2MS)

Williams, Kelleher, Heck

DOI: 10.1038/s41592-020-0764-5 DOI: 10.1021/acs.analchem.9b01669

Examples of Applications

Small molecules (metabolites/peptides)

Data averaging: sensitivity & quantitation accuracy

Full transients: ultra-high resolution, throughput

Transient post-processing: super-resolution algorithms

single ion counting





Large molecules (proteins)

Data averaging: top-down & middle-down sequencing

Transient post-processing: enhanced spectral quality

Transient post-processing:

Example: Isotopic Ratio Analysis of Uranium (UO₂)



Extreme spectral dynamic range is required (with high IR precision): 1: 20 000



I(²³⁸UO₂) 100 18 200





Simulation of FTMS data: DOI: 10.1021/jasms.0c00190

Averaging of the Unreduced Data: UO₂ ERSI



8000 scans data averaging

Markus et al., JASMS, in print







Challenge: quantitatively accurate analysis of diluted samples (low abundance)



Electron Impact



Increased sensitivity is provided by the unreduced data, improved LOD/LOQ



Examples of Applications

Small molecules (metabolites/peptides)

Data averaging: sensitivity & quantitation accuracy

Full transients: ultra-high resolution, throughput

Transient post-processing: super-resolution algorithms

single ion counting





Large molecules (proteins)

Data averaging: top-down & middle-down sequencing

Transient post-processing: enhanced spectral quality

Transient post-processing:



Upgraded Imaging Platform @ Leach Labs





LTQ Orbitrap XL (mFT) MALDI ion source Spectroglyph

Data acquisition system





Examples of Applications

Small molecules (metabolites/peptides)

Data averaging: sensitivity & quantitation accuracy

Full transients: ultra-high resolution, throughput

Transient post-processing: super-resolution algorithms

single ion counting





Large molecules (proteins)

Data averaging: top-down & middle-down sequencing

Transient post-processing: enhanced spectral quality

Transient post-processing:







Speed Up TMTc with Least-Squares Fitting (LSF)







Tsybin et al., Mass Spec. Rev. (2021) 10.1002/mas.21681

So, When to Use the Unreduced Data?



- Higher spectral dynamic range is required: data averaging across multiple scans
- Ultra-high resolution is required: ion detection during all the time ion signal rings
- Higher sensitivity is required: data averaging across multiple scans & multiple LC/GC-MS runs; data averaging of selected areas in single cell imaging Higher throughput is required: enhanced ion
- detection duty cycle due to parallel ion detection and ion manipulation (fragmentation, accumulation, etc.)
- Spectral complexity is beyond the conventional approaches: single ion counting, super-resolution

Always! But benefits differ between applications – to evaluate pros & cons



Spectroswiss Summary: Empowering FTMS with The Unreduced Data



Full profile

The two equal (!) types of the unreduced data in FTMS: transients and absorption mode mass spectra

eFT/mFT

- Advances in hardware and software enabled overcoming the technological challenges & creating the required tools for the unreduced data generation and handling
- Applied benefits: enhanced analytical characteristics for selected (especially multi-scan) applications & enabling the new ones: single ion counting, super-resolution MS, ...
- Current inhibitors: storage space for the unreduced data; lacksquarecomputational speed for super-resolution data processing

The unreduced data and tools are now available to advance FTMS applications

Full profile



Centroided

Empowering FTMS with The Unreduced Data Suggested Reading

- Resource on Orbitrap models design and applications: https://planetorbitrap.com/
- Advanced fundamentals in Fourier transform mass spectrometry: https://doi.org/10.1016/B978-0-12-814013-0.00005-3
- Enhanced Fourier transform for Orbitrap mass spectrometry: https://doi.org/10.1016/j.ijms.2014.05.019
- Absorption mode Fourier transform for FTMS: http://www.kilgourlab.com/absorption-mode-for-ft-ms/
- Transient-Mediated Simulations of FTMS Isotopic Distributions and Mass Spectra to Guide Experiment Design and Data Analysis: https://doi.org/10.1021/jasms.0c00190
- Fourier transform mass spectrometry at the uncertainty principle limit for improved qualitative and quantitative molecular analyses (PhD thesis, Anton Kozhinov): https://infoscience.epfl.ch/record/205045
- Trace-Level Persistent Organic Pollutant Analysis with Gas-Chromatography Orbitrap Mass Spectrometry Enhanced Performance by Complementary Acquisition and Processing of Time-Domain Data: https://doi.org/10.1021/jasms.9b00032
- Improved Uranium Isotope Ratio Analysis in the Liquid Sampling-Atmospheric Pressure Glow Discharge/Orbitrap FTMS Coupling Through the Use of an External Data Acquisition System: by Marcus et al., JASMS, in print

Questions? Ideas? Most welcome to discuss: tsybin@spectroswiss.ch



